



MESI

Simplifying Diagnostics



MESI ABPI MD

World's fastest Ankle-Brachial Index




What is Peripheral Arterial Disease (PAD)?

PAD is a circulatory problem where narrowed arteries reduce blood flow to the limbs.

The **global prevalence of PAD was 5.6%** in 2015, indicating that **≈236 million adults were living with PAD worldwide.**^[1]

The **principal cause of PAD is atherosclerosis** (>90% of cases).^{[2][3]} Since arteries in the legs are most commonly affected, the disease is also known as **LEAD** (Lower Extremity Artery Disease).



Risk factors


MANAGEABLE

- Smoking
- Diabetes mellitus
- Hypertension
- Hypercholesterolemia

NON-MANAGEABLE

- Age
- Genetic factors
- Gender
- Environment

Under-diagnosis of PAD in primary care may be a significant issue as **most patients with PAD do not experience typical symptoms.** The symptoms mainly appear during activities that increase energy demand.^[4]



Because of the high co-occurrence of PAD with Coronary Artery Disease (CAD) and Cerebrovascular Disease (CVD), **patients diagnosed with PAD have greater chance of early diagnosis of CAD and CVD.**

If not discovered and treated on time, PAD can lead to:



Heart attack



Stroke



Ulcerations




Amputations

How can PAD be detected?

Ankle-Brachial Index (ABI) is a reliable predictor of the occlusion of lower extremity arteries - PAD. Detection of PAD is even more important because **over 70% of the population is not aware that they have occluded arteries.**

ABI is a very simple **comparison of the blood pressures in the legs and the arms.** It is non-invasive and painless.



$$ABI = \frac{\text{SYSTOLIC BP}_{\text{ANKLE}}}{\text{SYSTOLIC BP}_{\text{ARM}}}$$

1.41 or more	1.40 - 1.00	0.99 - 0.91	0.90 - 0.51	0.50 or less
NON-COMPRESSIBLE	NORMAL	BORDERLINE	ABNORMAL	SEVERE
Indication of PAD/LEAD	Monitor risk factors/symptoms	Monitor risk factors/symptoms regularly	Indication of PAD/LEAD	Indication of PAD/LEAD
Perform advanced vascular testing or refer to a specialist			Perform advanced vascular testing or refer to a specialist	Perform advanced vascular testing or refer to a specialist

Measurement of Ankle-Brachial (pressure) Index or ABI offers **great accuracy and specificity in diagnosing PAD.** It is one of the most important tools for this purpose, especially due to its **non-invasive nature and cost-effectiveness.**

Areas of use

First-contact care and cardiovascular disease assessment

Every patient at risk for Peripheral Arterial Disease regardless of the presence of symptoms. ESC guidelines recommend early ABI measurement for the three groups listed below.

Patients with clinical suspicion for LEAD

- ▶ Unnoticeable pulse
- ▶ Claudication or symptoms suggestive of LEAD
- ▶ Non-healing wound

Patients with clinical conditions increasing LEAD risk

- ▶ CAD
- ▶ Heart failure
- ▶ Abdominal Aortic Aneurysm
- ▶ CKD

Asymptomatic individuals at risk for LEAD

- ▶ < 65 years with cardiovascular risk factors
 - Diabetes
 - Hypertension
 - Smoking
 - Dyslipidaemia
- ▶ < 50 years with family history of LEAD
- ▶ Everyone > 65 years


Compression therapy and wound care

Ankle-Brachial Index is a crucial measurement in wound care management and compression therapy.

ABI measurement should be performed to:

- ▶ Determine adequate arterial blood flow prior to compression therapy
- ▶ Rule out PAD/LEAD with a lower extremity wound
- ▶ Assess wound healing potential
- ▶ Evaluate therapeutic outcome
- ▶ Determine safe level of compression

RELATIONSHIP BETWEEN ABI AND COMPRESSION APPLIED ^[7]



Chronic disease management: Diabetes, Chronic kidney disease

Cardiovascular risk factors affect the onset of PAD differently, but always significantly. PAD is present in as many as **50% of patients** with diabetic foot ulceration (DFU).

DFU combined with PAD makes the wound **less likely to heal and more likely to require amputation** compared to patients without PAD.

Patients with impaired renal function have **greater than two-fold risk of developing PAD**. PAD is usually not prognosed well and often coexists with other conditions like CAD and diabetes. In end-stage renal disease, **diabetes is the most common factor** for PAD risk and outcomes.


MESI ABPI MD

The MESI ABPI MD device is an **automated Ankle-Brachial Index system**. With its **3CUFF™ technology**, it allows simultaneous measurement of brachial and ankle pressures. It also incorporates the **PADsense™ algorithm** for the detection of severe Peripheral Arterial Disease.



Highlights of MESI ABPI MD:


- ▶ Guided measurement procedure on the device screen
- ▶ SmartArm™ detection – to determine the higher blood pressure of the two
- ▶ Automated, 3-cuff simultaneous measurement
- ▶ Multiple cuff sizes and ability to mix-and-match different size cuffs during one measurement
- ▶ Discovers blocked arteries in 1 minute
- ▶ Pulse waveforms and oscillation graphs
- ▶ Advanced review and alerts thanks to PADsense™ algorithm



Pulse waveforms and oscillation graphs

The MESI ABPI MD uses the PADsense™ pattern recognition algorithm to **automatically interpret the acquired pulse waveform** and calculate the ABI with the result.

Combining the ABI result and pulse waveform represents **the best practice in evaluating the presence and severity** of Peripheral Arterial Disease (PAD).



NORMAL RESULT

- 1 A rapid rise in the upstroke during systole
- 2 A very sharp peak
- 3 A gradual downstroke
- 4 A presence of dichototic notch


ABNORMAL RESULT

A flattened pulsewave recording or one without the typical lemon shape is an indicator of severe PAD.

The absence of the pulsations caused by occlusions in the artery makes it impossible to calculate the ankle pressures. Instead of the ABI value, the device will display a “PAD” result.

The measurement of the Ankle-Brachial Index also provides **oscillation graphs, which provide the analysis of the entire deflation portion** of the measurement.

The oscillation graph in a healthy patient has its **easily recognizable lemon shape, which shows the arteries are elastic**. If the patient has PAD or incompressible arteries, the graph will reflect the condition.



Normal oscillation graph



Severe PAD oscillation graph



Incompressible arteries oscillation graph

Technical specifications

Measurements

Ankle-Brachial Pressure Index using improved oscillometric method and plethysmography, with PADsense™ algorithm

Heart rate and systolic, diastolic, and mean blood pressure using an improved oscillometric method and plethysmography

Measurement extensions

Blood pressure

Pulse wave velocity*

*Available with MESIresults software

Device dimensions

Width: 223 mm / 8.78 in

Depth: 174 mm / 6.85 in

Height: 86 mm / 3.38 in

Weight: 1000 g / 2.2 lbs

Power & battery

Battery type: Rechargeable lithium-polymer battery

Capacity: 4400 mAh

AC/DC adaptor: FRIWO FW8030M/05 or FOX NEO30-XM

Input: 100-240 V AC / 50-60 Hz / 600-300 mA

Output: 5V DC / 3.0 A

Examinations per battery charge: > 100

Charge time for depleted battery: 2 hrs

Touchscreen display

Diagonal size: 8"

Resolution: 800 x 600 pixels

Color depth: 16-bpp

Cuffs

Medium size cuffs (circumference: 22-32 cm / 6.7-12.6 in)

Large size cuffs: (circumference: 32-42 cm / 12.6-16.5 in)

Measurement range & Accuracy

Measurement range:

Pressure: 0 to 299 mmHg

Pulse rate: 30 to 199 beats per minute

Accuracy:

Pressure: ± 3 mmHg

Heart rate: $\pm 5\%$ of value

Ankle-Brachial Pressure Index: ± 0.1

Data backup and management

MESIresults software

Protection Classification

Protection against electric shock: Class II

Medical device classification: Class IIa

Applied parts (cuffs for arms and ankles):

Type BF Applied part

Software classification: Class A

RF emissions (CIPSR 11): Group 1, Class A

Ingress protection: IP2X

Applied Standards

EN 60601-1:2006+A1:2013 General requirements for basic safety and essential performance

EN 60601-1-2:2015+A1:2021 Electromagnetic disturbances - Requirements and tests

EN 80601-2-30:2019 Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers

Operating Conditions

Temperature, operating: 10° to 40°C (50° to 104°F)

Relative humidity: 30 to 80% (no condensation)

Pressure during operation: 700 to 1060 hPa

Transport & Storage Conditions

Temperature: 0° to 40°C (32° to 104°F)

Relative humidity: 15 to 85%

Pressure during storage: 700 to 1060 hPa

Measure ABI with MESI ABPI MD in three simple steps

Place the colour-coded cuffs
on the arms and ankles **1**



2 Press the **START** button for
the measurement to begin



Analyse the pulse waveforms
and read the ABI results **3**

References

1. Song P, Fang Z, Wang H, Cai Y, Rahimi K, Zhu Y, Fowkes FGR, Fowkes FJI, Rudan I. Global and regional prevalence, burden, and risk factors for carotid atherosclerosis: a systematic review, meta-analysis, and modelling study. *Lancet Glob Health*. 2020; 8:e721–e729. doi: 10.1016/S2214-109X(20)30117-0
2. Hoyer C, Sandermann J, Petersen LJ. The toe-brachial index in the diagnosis of peripheral arterial disease. *J Vasc Surg*. 2013;58:231–238
3. Thukkani AK, Kinlay S. Endovascular intervention for peripheral artery disease. *Circ Res*. 2015;116:1599–1613
4. Jelani QU, Petrov M, Martinez SC, Holmvang L, Al-Shaibi K, Alasnag M. Peripheral Arterial Disease in Women: an Overview of Risk Factor Profile, Clinical Features, and Outcomes. *Curr Atheroscler Rep*. 2018 Jun 02;20(8):40
5. Hirsch AT, Criqui MH, Treat-Jacobson D, et al. Peripheral Arterial Disease Detection, Awareness, and Treatment in Primary Care. *JAMA*. 2001;286(11):1317–1324. doi: 10.1001/jama.286.11.1317
6. McDermott MM, Greenland P, Liu K, et al. Leg Symptoms in Peripheral Arterial Disease: Associated Clinical Characteristics and Functional Impairment. *JAMA*. 2001;286(13):1599–1606. doi: 10.1001/jama.286.13.1599
7. Vowden P and K (2001). Doppler assessment and ABPI: Interpretation in the management of leg ulceration. *World Wide Wounds*. <http://www.worldwidewounds.com/2001/march/Vowden/Doppler-assessment-and-ABPI>

MESI, Ltd.
Leskoškova cesta 11a
1000 Ljubljana, Slovenia, EU

 MESIdoo

 MESImedical

 MESI



European production
and development



Worldwide presence

FDA

Food and Drug
Administration
cleared



EU Medical Device
Regulations
compliant



ISO 9001 Q-1664
ISO 13485 M-049

ISO 9001 and
ISO 13485 certified



MDSAP compliant